DE: Vor Sonnenlicht geschützt aufbewahren
EN: Protect from sunlight
FR: Ne pas exposer aux rayons solaires
ES: Mantener fuera de la luz solar
IT: Non esporre a raggi solari
PT: Guardar em lugar seco e protegido do sol
DA: Opbevares væk fra sollys / SV: Utsätt inte för direkt solljus
NO: Oppbevares unna sollys / FI: Säilytettävä auringonvaloista suojattuna
TR: Güneş ışıklarından uzak saklayın
EL: Αποθήκευτο χωρίς έκθεση στο ηλιακό φώς
UK: Зберігати в захищенному від сонячних променів місці
RU: Не допускать попадання сонячних лучей
AR: يحفظ بعيداً عن الشمس


DE: Trocken aufbewahren
EN: Store in dry place
FR: Conserver au sec / ES: Mantener seco
IT: Teme l'umidità / PT Guardar em lugar seco
DA: Opbevares tørt / SV: Förvaras tørt / NO: Oppbevares tørt
FI: Säilytettävä kuivassa
TR: Kuruyerde saklayın / EL: Να φυλασσοται σε ξηρό μέρος
UK: Зберігати сухим місцем / RU: Хранить в сухом месте
AR: تحفظ في مكان جاف

DE: Hersteller
EN: Manufacturer
FR: Fabricant / ES: Fabricante
IT: Fabricante / PT: Fabricante
DA: Produsent / SV: Tillverkare / NO: Produsent / FI: Valmistaja
TR: Üretici / EL: Κατασκευαστής
UK: Виробник / RU: Изготовитель
AR: مصنّع

DE: Konformitätszeichen nach EU Richtlinie 93/42/EWG
EN: Conformity mark according to EU Council Directive 93/42/EEC
FR: Marque de conformité selon la directive communautaire 93/42/CEE
ES: Marca de conformidad según la directiva comunitaria 93/42/CEE
IT: Marchio CE secondo la direttiva 93/42/CEE
PT: Marcção de conformidade de acordo com a Diretiva da UE 93/42/CEE
DA: Overensstemmelsesmærke i henhold til EU-direktiv 93/42/EØF
SV: Översensstämmelsemärke enligt EU-direktiv 93/42/EEG
NO: Samensvermærke i henhold til EU-direktiv 93/42/EEC
FI: Vaatimustenmukaisuusmerkintä EY-direktiivin 93/42/ETY mukaisesti
TR: 93/42 EWG sayılı AB direktifin gereğince uygunluk beyanı
EL: Σήμα καταληξτάς σύμφωνα με την οδηγία της ΕΕ 93/42/EOK
UK: Знак відповідності відповідно до Директиви ЄС 93/42/ЕС
RU: Знак соответствия согласно Директиве ЕС 93/42/EС
AR: علامات التوافق طبقاً لـ إرشادات الاتحاد الأوروبي 93/42/EC الجمدة الأوروبية


DE: Gebrauchsanweisung befolgen
EN: Observe operating instructions
FR: Respectez le mode d'emploi
ES: Siga las instrucciones de empleo
IT: Rispettare le istruzioni per l'uso
PT: Favor de observar o manual de operação
DA: Følg brugsanvisning / SV: Följ bruksanvisning
NO: Følg bruksanvisning / FI: Noudata käyttöohjeita
TR: Kullanma kilavuzunu dikkate alınır
EL: Ακολουθήστε τις οδηγίες χρήσης
UK: Дотримуватися інструкції з експлуатації
RU: Следуйте указаниям в инструкции по эксплуатации
AR: يجب اتباع تعليمات الاستخدام


DE: Bestellnummer
EN: Catalog number
FR: Référence du catalogue / ES: Número de catálogo
IT: Codice prodotto / PT: Número de enconda
DA: Löbenummer / SV: Beställningsnummer
NO: Ordrenummer / FI: Tilauksennumero
TR: Sipariş numarası / EL: Αριθμός παραγγελίας
UR: Номер для замовлення / RU: Номер для заказа
AR: رقم المطلب


DE: Chargenbezeichnung
EN: Batch code
FR: Code du lot / ES: Código de lote
IT: Numero del lotto / PT: Designação do lote
DA: Batchetegnelse af partiet / SV: Batchteckning
NO: Batchtegnesle / FI: Erän nimittys
TR: Parti numarası / EL: Στοχεύο τοποτοποίησης της παρτίδας
UK: Номер партії / RU: Обозначение партии Номер партії
AR: رقم التشغيل


Small D = 18 mm
REF 281481
PZN 17373455

Medium D = 22 mm
REF 281482
PZN 17373461


Small D = 18 mm
REF 281441
PZN 17373389

Medium D = 20 mm
REF 281448
PZN 17373403


Large D = 28 mm
REF 281442
PZN 17373426

1.

2.

3.

4.

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1. Purpose**Why using mamivac® nipple shields?**

Many mothers use conventional round shaped nipple shields when having sensitive and sore nipples. These should be used for a short time only, because they have influence on the natural way of breastfeeding. The round shape covers mum's skin, and it may happen that the baby refuses to suck when mum is trying to breastfeed again without a nipple shield. This is because the baby does not identify its mother due to the irritation by the nipple shield. The special shape of **mamivac®** nipple shields avoids this by reason of not interrupting the contact between mother and baby. Therefore the baby reminds the natural smell of its mother.

2. Use

Clean your hands before using **mamivac®** nipple shields. A short massage has the effect that nipples stand upright and make the use of **mamivac®** nipple shields much easier. Turn the soft edges up and press the **mamivac®** nipple shield as close as possible to your nipple. Press the soft edges afterwards to your skin. Wetting the **mamivac®** nipple shield supports the attachment to your skin. Make sure to place the **mamivac®** nipple shield with the flat side up that your baby's nose has direct contact to your skin while breastfeeding. Then the baby identifies your smell and feels safe. Do not use Vaseline or other creams. For further information regarding the use please ask your lactation consultant or midwife. Prior to using the nipple shields, check them for cracks and proper functionality.

3. Cleaning and disinfection

mamivac® nipple shields consist of thin, transparent and tasteless silicone. Prior to the first use and before any other use, disinfect the nipple shields by boiling them in water for at least 5 minutes. Clean the nipple shields after each use with commercial washing-up liquid and rinse them thoroughly with lukewarm water. With proper use and storage, the **mamivac®** nipple shields have a service life of 2 years. Keep the operating instructions. This product is not a toy. Keep it away from children.

4. Disposal

Please dispose of your **mamivac®** nipple shields in accordance with your local disposal regulations for plastic products.

</