

FEEDING SUPPORT

mamivac® BABYFEED

- Practical infant feeding helper for babies
- Facilitates drinking for babies
- Does not stimulate the baby's sucking instinct
- Can be disinfected by boiling water
- Reusable
- Made of polypropylene



Content: 2 pieces **REF 282112**
Clinic pack: 10 pieces **REF 282113**

mamivac® BABYCUP

Safe and controlled administration of precious colostrum

- A unique, simple tool for the administration of collected breast milk
- Allows feeding at the baby's natural pace
- No leaking or spilling
- Also for cleft lip and/or palate (cheiloschisis)
- Fits all bottles with standard thread
- With adapter, it can also be used for wide-neck bottles



Babycup
REF 281610

Bottle 50 ml
REF 2801135

Babycup-Set
with wide-neck adapter
and 140 ml bottle
REF 281611



SORE NIPPLES

mamivac® NIPPLE SHIELDS

- To be used on sensitive and sore nipples
- Ergonomic design
- Best body contact between mother and child
- Can be disinfected by boiling in water
- Made of thin, transparent and tasteless silicone
- Help also with flat and retracted nipples or with baby's nipple confusion

»»cherry««

Small (18 mm)
REF 281481

Medium (22 mm)
REF 281482

Large (28 mm)
REF 281515

»conical«

Small (18 mm)
REF 281441

Medium (20 mm)
REF 281448

Large (28 mm)
REF 281442

Clinic pack with 50 pcs. available

mamivac® BREAST COMPRESS

- Individually sterile packed dressing pad for the treatment of sore nipples
- Pleasant wearing comfort
- Absorbs wound secretion
- Supports wound healing without active ingredients
- Low-adherent wound compress



1 Unit = 10 pieces
REF 281460

mamila® BREAST CREME

- Sugar and paraffin free
- Free Of preservatives
- Nurtures rough nipples
- Does not have to be washed off prior to breastfeeding
- Contains pure lanolin
- Ideal for sensitive and dry skin
- Can prevent pregnancy stretch marks or can also be used as lip balm



Tube 30 ml **REF 281488**

FOR THE BREAST

mamivac® NURSING PADS

- The antimicrobial nursing pad reduces the growth of germs like bacteria, fungi and yeast.
- Antimicrobial, highly absorbent core
- Non-leaking, clothes remain dry
- The antimicrobial effect is not washed out.
- Breathable microfibre keeps skin dry.
- Especially skin-friendly, soft and safe
- No unpleasant odours
- Welded and tight, yet soft seams
- Anatomically shaped for best wearing comfort
- Washable up to 60 °C, reusable

6 pieces
with washing bag
in high-grade box
REF 281495



antimicrobial
(silver ions)



highly absorbent
and non-leaking




breathable

mamivac® THERMO-PACK

The **mamivac® Thermo-Pack** is an anatomically shaped gel cushion which can be used to either warm or cool the breast.

- Includes 3 hygienic covers
- Assists with many breastfeeding problems
- For optimal relaxation
- Universal use



1 piece with
3 hygienic covers
REF 281487



SPECIAL EQUIPMENT

CBC FUNwatch

Nurse watch with convenient carabiner.



REF 6098 REF 6093 REF 6096 REF 6091 REF 6008 REF 6187

CBC CLIPwatch



REF 6177 REF 6178 REF 6186

Domings for CLIPwatch



Puls Red Blue Silver Heart Mother

mamivac® PINARD STETHOSCOPE

- Turned beech wood
- Available in two lengths

Pinard short 18 cm **REF 2260**
Pinard long 33 cm **REF 2261**



€ Made in Germany

mamivac® FOR MOTHERS. FOR BABIES.



BREASTFEEDING PRODUCTS

mamivac.com
Made in Germany

BREASTFEEDING –
most beautiful time
for mother and baby

Breast milk is the best food you can
give your baby because it contains
everything your baby needs to prosper.

However, there may be situations
that can affect the breastfeeding
relationship between mother
and child.

As a manufacturer of
medical devices, we at
Kimetec are certified
in accordance with
ISO 13485 and would
like to support you with
high-quality, safe and
effective products.

We wish you an enjoyable and happy breastfeeding time.

WHICH PUMP FOR WHAT PURPOSE?

We have the right pump for you. Select the pump which is most
appropriate to meet your personal expectations.

Medical reasons:	mastitis, sore nipples, flat or retracted nipples, plugged ducts, insufficient lactation, baby in hospital
Use in hospitals and rentals:	full-time occupation, baby in hospital, premature infants, frequent use
Private use:	full-time occupation, frequent use, baby in hospital
Occasional use:	on the go, part-time occupation, strain in the breast, shifts in breastfeeding times

Manual breast pump mamivac® EASY	
Electric breast pump mamivac® LACTIVE	
Electric breast pump mamivac® SENSITIVE-C	

TECHNOLOGY

mamivac®
2-PHASE TECHNOLOGY

To make the expressing of breast milk as
natural as possible, we have developed
the mamivac® 2-phase technology.
The continuously adjustable phases
consist of a stimulation phase and an
expressing phase.

The stimulation phase imitates
the fast and weak sucking of the
baby at the beginning and gently
massages the breast to stimulate
the flow of milk. The gentle and effective
expressing phase simulates the slower,
more intense sucking of the baby to
efficiently express the invaluable breast milk.



mamivac®
DOUBLE PUMPING

The simultaneous expressing of milk
from both breasts saves time
reducing the expressing time by about 50 %.
This gives the mom more time – not
only with multiple births. The patented
valve technology ensures double
pumping without suction loss.
All our electric breast pumps offer
direct connections for two sets.

mamivac®
PROTECTIVE DIAPHRAGM

The patented protective diaphragm in
the suction kit offers the best possible
protection and safety for mother
and baby.


This diaphragm hermetically seals off
the breast pump creating a closed
system. Breast milk, bacteria or
germs cannot contaminate neither
the hoses nor the pump. This avoids
cross contamination between
users and the device itself.



ELECTRICAL BREAST PUMPS

mamivac®
SENSITIVE-C For rental use.

- The mamivac® 2-phase technology.
imitates the baby's sucking
- Continuously
adjustable
vacuum strength,
freely selectable
pumping interval
- Vacuum strength
corresponds
to pumps used
in hospitals
- No suction loss with
double pumping
- Gentle, hygienic
and silent operation guaranteed
- Compact, maintenance-free
and durable device
- Germ free – patented protective
diaphragm in the suction kit offers the
best safety for mother and child



mamivac®
LACTIVE
The compact 2-phase breast pump
for active and working mothers.

- Easy to use
- Expressing comfort due to
2-phase technology
- Visual stimulation with a photo
of your baby. This stimulus
supports the flow of breast milk.
- Time-saving due to optional
double pumping
- Best hygienic standard
for mother & child



REF 287000
Made in Germany

ACCESSORY

mamivac®
CARRYING CASE

Robust case for safe transportation
and storage of the electrical breast
pump mamivac® SENSITIVE-C.

REF 283011



mamivac®
COOL BAG

The expressed milk can be
transported easily and safely
in this cooling bag.
The spacious compartment can
hold the mamivac® EASY and
LACTIVE and additional bottles.

REF 283024



mamivac®
SUCTION KIT PREMIUM

- Free of bisphenol A with wide neck
bottle for maximum comfort
- Easy and safe operation
- Protective diaphragm ensures best
possible safety for mother and child
- Fits the breast pumps of the mamivac®
SENSITIVE and LACTIVE series and
other piston breast pumps
- Silicone soft insert (Diameter 24 mm)
only for SUCTION KIT PREMIUM L

Breast shield in different sizes:
S Diameter 20 mm
L Diameter 28 mm
XL Diameter 36 mm



SUCTION KIT PREMIUM S
SUCTION KIT PREMIUM L Single Kit
SUCTION KIT PREMIUM L Double Kit
SUCTION KIT PREMIUM XL

REF 283129
REF 283120
REF 283123
REF 283133

MANUAL BREAST PUMP

mamivac®
EASY The simple solution when on the go.

- As powerful as an electric breast pump
- The mamivac® 2-phase technology
imitates the baby's sucking.
- Adjustment valve for continuously
adjustable suction strength
- Consists of a few parts only
– easy cleaning
- Comfortable and safe
one-hand operation
- Soft silicone breast shield
- Hygienic as the invaluable
breast milk flows straight
into the bottle



REF 281360
Made in Germany



mamivac®
BREAST MILK FREEZER BAGS

For a hygienic freezing of your invaluable
breast milk.

- Integrated closure
- Safely closed – the milk can be poured
out again without spilling
- With scaling, can be labelled
- Safe stand due to practical
folding bottom



1 Unit = 20pieces
REF 282100

Kimetec GmbH
Gerlinger Strasse 36-38
71254 Ditzingen, Germany

Tel.: +49 (0) 71 56 / 1 76 02 - 200
Fax.: +49 (0) 71 56 / 1 76 02 - 500
E-mail: info@kimetec.de

mamivac®.com
Made in Germany