

FEEDING SUPPORT

mamivac® BABYFEED

- Practical infant feeding helper for babies
- Facilitates drinking for babies
- Does not stimulate the baby's sucking instinct
- Can be disinfected by boiling water
- Reusable
- Made of polypropylene



Content: 2 pieces **REF 282112**
Clinic pack: 10 pieces **REF 282113**



mamivac® BABY CUP

Safe and controlled administration of precious colostrum

- A unique, simple tool for the administration of collected breast milk
- Allows feeding at the baby's natural pace
- No leaking or spilling
- Also for cleft lip and/or palate (cheiloschisis)
- Fits all bottles with standard thread
- With adapter, it can also be used for wide-neck bottles



Babycup
REF 281610

Bottle 50 ml
REF 2801135

Babycup-Set
with wide-neck adapter
and 140 ml bottle
REF 281611

SORE NIPPLES

mamivac® NIPPLE SHIELDS

- To be used on sensitive and sore nipples
- Ergonomic design
- Best body contact between mother and child
- Can be disinfected by boiling in water
- Made of thin, transparent and tasteless silicone
- Help also with flat and retracted nipples or with baby's nipple confusion

»»cherry««
Small (18 mm)
REF 281481

Medium (22 mm)
REF 281482


Large (28 mm)
REF 281515



»conical«
Small (18 mm)
REF 281441

Medium (20 mm)
REF 281448

Large (28 mm)
REF 281442



Clinic pack with 50 pcs. available

mamivac® BREAST COMPRESS

- Individually sterile packed dressing pad for the treatment of sore nipples
- Pleasant wearing comfort
- Absorbs wound secretion
- Supports wound healing without active ingredients
- Low-adherent wound compress



1 Unit = 10 pieces
REF 281460

mamivac® BREAST CREME

- Sugar and paraffin free
- Free of preservatives
- Nurtures rough nipples
- Does not have to be washed off prior to breastfeeding
- Contains pure lanolin
- Ideal for sensitive and dry skin
- Can prevent pregnancy stretch marks or can also be used as lip balm



Tube 30 ml **REF 281488**

FOR THE BREAST

mamivac® NURSING PADS

- The antimicrobial nursing pad reduces the growth of germs like bacteria, fungi and yeast.
- Antimicrobial, highly absorbent core
- Non-leaking, clothes remain dry
- The antimicrobial effect is not washed out.
- Breathable microfibre keeps skin dry.
- Especially skin-friendly, soft and safe
- No unpleasant odours
- Welded and tight, yet soft seams
- Anatomically shaped for best wearing comfort
- Washable up to 60 °C, reusable

6 pieces
with washing bag
in high-grade box
REF 281495



antimicrobial
(silver ions)



highly absorbent
and non-leaking



breathable

mamivac® THERMO-PACK

The **mamivac® Thermo-Pack** is an anatomically shaped gel cushion which can be used to either warm or cool the breast.

- Includes 3 hygienic covers
- Assists with many breastfeeding problems
- For optimal relaxation
- Universal use



1 piece with
3 hygienic covers
REF 281487



SPECIAL EQUIPMENT

CBC FUNWATCH

Nurse watch with convenient carabiner.



REF 6098 REF 6093 REF 6096 REF 6091 REF 6008

CBC CLIPWATCH



REF 6177 REF 6178 REF 6186

Domings for ClipWatch



Puls



Heart



Empathy



Mother

mamivac® PINARD STETHOSCOPE

- Turned beech wood
- Available in two lengths

Pinard short 18 cm **REF 2260**
Pinard long 33 cm **REF 2261**



CE Made in Germany

mamivac®
FOR MOTHERS. FOR BABIES.



BREASTFEEDING PRODUCTS

mamivac.com
Made in Germany

BREASTFEEDING – most beautiful time for mother and baby

Breast milk is the best food you can give your baby because it contains everything your baby needs to prosper.

However, there may be situations that can affect the breastfeeding relationship between mother and child.

As a manufacturer of medical devices, we at Kimetec are certified in accordance with ISO 13485 and would like to support you with high-quality, safe and effective products.

We wish you an enjoyable and happy breastfeeding time.

WHICH PUMP FOR WHAT PURPOSE?

We have the right pump for you. Select the pump which is most appropriate to meet your personal expectations.

Medical reasons:	mastitis, sore nipples, flat or retracted nipples, plugged ducts, insufficient lactation, baby in hospital
Use in hospitals and rentals:	full-time occupation, baby in hospital, premature infants, frequent use
Private use:	full-time occupation, frequent use, baby in hospital
Occasional use:	on the go, part-time occupation, strain in the breast, shifts in breastfeeding times

Manual breast pump **mamivac® EASY**

Electric breast pump **mamivac® LACTIVE**

Electric breast pump **mamivac® SENSITIVE-C**



TECHNOLOGY

mamivac® 2-PHASE TECHNOLOGY


To make the expressing of breast milk as natural as possible, we have developed the **mamivac® 2-phase technology**. The continuously adjustable phases consist of a stimulation phase and an expressing phase.

The **stimulation phase** imitates the fast and weak sucking of the baby at the beginning and gently massages the breast to stimulate the flow of milk. The gentle and effective **expressing phase** simulates the slower, more intense sucking of the baby to efficiently express the invaluable breast milk.



mamivac® DOUBLE PUMPING


The simultaneous expressing of milk from both breasts **saves time** reducing the expressing time by about 50%. This gives the mom more time – not only with multiple births. The patented valve technology ensures double pumping without suction loss. All our electric breast pumps offer direct connections for two sets.



mamivac® PROTECTIVE DIAPHRAGM

The patented protective diaphragm in the suction kit offers the **best possible protection and safety for mother and baby**.


This diaphragm hermetically seals off the breast pump creating a closed system. Breast milk, bacteria or germs cannot contaminate neither the hoses nor the pump. This avoids cross contamination between users and the device itself.



ELECTRICAL BREAST PUMPS

mamivac® SENSITIVE-C For rental use.

- The **mamivac® 2-phase technology** imitates the baby's sucking
- Continuously adjustable vacuum strength, freely selectable pumping interval
- Vacuum strength corresponds to pumps used in hospitals
- No suction loss with double pumping
- Gentle, hygienic and silent operation guaranteed
- Compact, maintenance-free and durable device
- Germ free – patented protective diaphragm in the suction kit offers the best safety for mother and child



REF 284000
Made in Germany

mamivac® LACTIVE The compact 2-phase breast pump for active and working mothers.

- Easy to use
- Expressing comfort due to 2-phase technology
- Visual stimulation with a photo of your baby. This stimulus supports the flow of breast milk.
- Time-saving due to optional double pumping
- Best hygienic standard for mother & child



REF 287000
Made in Germany

ACCESSORY

mamivac® CARRYING CASE

Robust case for safe transportation and storage of the electrical breast pump **mamivac® SENSITIVE-C**.

REF 283011



mamivac® COOL BAG

The expressed milk can be transported easily and safely in this cooling bag. The spacious compartment can hold the **mamivac® EASY** and **LACTIVE** and additional bottles.

REF 283024



mamivac® SUCTION KIT PREMIUM

- Free of bisphenol A with wide neck bottle for maximum comfort
- Easy and safe operation
- Protective diaphragm ensures best possible safety for mother and child
- Fits the breast pumps of the **mamivac® SENSITIVE** and **LACTIVE** series and other piston breast pumps
- Silicone soft insert (Diameter 24 mm) only for SUCTION KIT PREMIUM L

Breast shield in different sizes:
S Diameter 20 mm
L Diameter 28 mm
XL Diameter 36 mm

SUCTION KIT PREMIUM S
SUCTION KIT PREMIUM L Single Kit
SUCTION KIT PREMIUM L Double Kit
SUCTION KIT PREMIUM XL



REF 283129
REF 283120
REF 283123
REF 283133

MANUAL BREAST PUMP

mamivac® EASY The simple solution when on the go.

- As powerful as an electric breast pump
- The **mamivac® 2-phase technology** imitates the baby's sucking.
- Adjustment valve for continuously adjustable suction strength
- Consists of a few parts only – easy cleaning
- Comfortable and safe one-hand operation
- Soft silicone breast shield
- Hygienic as the invaluable breast milk flows straight into the bottle



REF 281360
Made in Germany

mamivac® BREAST MILK FREEZER BAGS

For a hygienic freezing of your invaluable breast milk.

- Integrated closure
- Safely closed – the milk can be poured out again without spilling
- With scaling, can be labelled
- Safe stand due to practical folding bottom



1 Unit = 20pieces
REF 282100

Kimetec GmbH
Dornierstr. 27
71254 Ditzingen, Germany

Tel.: +49 (0) 71 56 / 1 76 02 - 200
Fax.: +49 (0) 71 56 / 1 76 02 - 500
E-mail: info@kimetec.de

mamivac®.com
Made in Germany