


# mamivac<sup>®</sup>

 MAMILLOMETER

## Find the correct size of nipple shields




conical

S = 18 mm

M = 20 mm






L = 28 mm



cherry shaped

S = 18 mm



M = 22 mm



# How do you find the correct size?

Place the different punched rings over your nipple and check which size you need.

- The size of the nipple shield depends on the mamilla and not on the areola.
- Take the measurement in a relaxed atmosphere.  
If the nipple has about 2 mm of leeway, you have found the nipple shield with the correct fit.
- If you are not sure, ask your lactation consultant for support.



## Our breast shield and suction kit sizes

